



**PROCEDURE FOR THE DETERMINATION OF POLYNUCLEAR
AROMATIC HYDROCARBONS IN VEHICULAR EXHAUST USING GAS
CHROMATOGRAPHY/MASS SPECTROMETRY**

**SOP MV-AEROSOL-144 Version 2.2
Effective Date: September 1, 2019**

Aerosol Analysis and Methods Evaluation Section
Chemical Analysis and Emissions Research Branch
Mobile Source Laboratory Division

DISCLAIMER: This procedure has been reviewed by the staff of the California Air Resources Board and approved for publication. Mention of any trade name or commercial product in this Standard Operating Procedure does not constitute endorsement or recommendation of this product by Air Resources Board. Specific brand names and instrument descriptions listed in the Standard Operating Procedure are for equipment used by the Air Resources Board laboratory.

Table of Contents

1.	Scope.....	1
2.	Summary of Test Method.....	1
3.	Terminology (Definitions).....	1
4.	Limitations and Interferences.....	2
5.	Safety.....	2
6.	Equipments and Apparatus.....	3
7.	Reagents and Materials.....	4
8.	Preparation of Apparatus	5
9.	Preparation of sampling materials.....	5
10.	Preparation of Standards for PAH Analysis.....	6
11.	Procedures for PAH Quantification.....	9
12.	Quality Control (QC)	12
13.	Calculation.....	15
14.	Report	17
15.	Reference Documents	17
16.	Revision History	17

Standard Operating Procedure for the Determination of Polynuclear Aromatic Hydrocarbons in Vehicular Exhaust using Gas Chromatography/Mass Spectrometry

1. Scope

- 1.1 This Standard Operating Procedure (SOP) describes the identification and quantification of polynuclear aromatic hydrocarbons (PAHs) in motor vehicle exhaust. The PAHs are extracted from vehicle exhaust collected on filters or with XAD/polyurethane foam (PUF) and are analyzed using gas chromatography/mass spectrometry (GC/MS).
- 1.2 The target analytes (PAHs compounds) in this GC/MS method are listed in Table 1.
- 1.3 The individual PAH detection limit for the extract is equal to or less than 0.1 ng/mL depending on the extent of interference and the sensitivity of each PAH's. The detection limit for the mass of each PAH on the filter and in the PUF/XAD/PUF cartridge can range from 0.12 ng to 10 ng depending upon the solvent purity, background interferences, and analytical conditions (variation of final extraction volume and GC injection volume).

2. Summary of Test Method

- 2.1 PAHs from vehicular emissions are collected on different types of sample media, including Teflon coated glass fiber filter, quartz fiber filter, polyurethane foam, and XAD resins.
- 2.2 PAHs on the sample media are solvent extracted and analyzed by gas chromatography/mass spectrometry (GC/MS).
- 2.3 Prior to extraction, known amounts of recovery standards (RS) are added to the sample. The sample is then extracted either with an accelerated solvent extraction (ASE) or a Soxhlet extractor using dichloromethane.
- 2.4 PAHs are identified by their retention times on GC and characteristic mass peaks. The recovery standards and target analytes are quantified using the calibration curves of each PAHs with corrections by response of internal standard.
- 2.5 The PAH mass (ng) is calculated by multiplying the extract volume and its concentration, and is corrected for its recovery efficiency.
- 2.6 A flow chart of this method is shown in Figure 1.

3. Terminology (Definitions)

- 3.1 Amu (Atomic mass unit) - Unit of mass to express atomic or molecular masses.
- 3.2 Internal Standard (IS) – Internal Standard is usually a deuterated compound added to a sample extract in a known amount and is used to calibrate concentrations of other analytes. The internal standard must not be the target analytes.

- 3.3 Method Detection Limit (MDL) – the minimum concentration of a substance that can be measured and reported with confidence and the value is above zero.
- 3.4 Polynuclear Aromatic Hydrocarbons (PAHs) – aromatic hydrocarbons with two or more fused aromatic rings.
- 3.5 Recovery Standard (RS) – Recovery Standard is usually a deuterated compound added to a sample in known amount prior to solvent extraction. The measured amount of the recovery standard after extraction is compared to its theoretical (calculated) value to establish the extraction efficiency and a correction factor for a specified analyte.

4. Limitations and Interferences

- 4.1 The target analytes are not limited to those listed in Table 1. Other compounds may also be added, depending on the availability of National Institute of Standards and Technology (NIST) reference materials or traceable standards.
- 4.2 The GC/MS analytical method is not limited to quantify PAHs collected on filter or in PUF/XAD/PUF cartridge. With proper extraction and sample preparation, the GC/MS method can also identify and quantify PAHs from other sample collection media and emission sources.
- 4.3 Quantification of the target analyte is subject to interference from compounds that have the same m/z with the same retention time region in the analysis.
- 4.4 Closely eluting compounds with the same quantification ions may not have enough peak resolution for quantification.
- 4.5 Teflon filters can only be extracted by Soxhlet extraction. The Accelerated Solvent Extraction (ASE) dissolves the support ring and causes contamination.

5. Safety

- 5.1 Many chemicals, especially PAHs, are toxic and/or carcinogenic. These chemicals must be handled extremely carefully with proper protection in the fume hood.
- 5.2 Solvents such as acetone and hexane are flammable and harmful. Handle these solvents in the fume hood.
- 5.3 Dichloromethane (DCM) is a potential carcinogen and can penetrate gloves. Handle DCM in the fume hood and avoid direct contact with skin.
- 5.4 The ASE is under high pressure (~ 1,600 psi) during operation and must be operated in the fume hood.

6. Equipments and Apparatus

6.1 Accelerated Solvent Extractor (or pressurized fluid extractor)

The PAH extraction is performed with the Accelerated Solvent Extractor (ASE, Dionex Model 350). The extractor is capable of using different sizes of extraction cells and allows automatic mixing and delivery of up to two separate solvents.

6.2 Soxhlet Extraction

The Soxhlet extraction apparatus is shown in Figure 2. The extractor must include the following components:

6.2.1 A chiller to control the coolant temperature to as low as 5 °C.

6.2.2 A heating mantle with a Variac to control solvent distillation rate.

6.2.3 A magnetic stirring bar or boiling chips to prevent solvent from superheating.

6.3 RapidVap N₂ Dry Evaporation System

The RapidVap N₂ Dry Evaporation System use a combination of nitrogen flow down, dry heat and vortex motion to evaporate solvent of liquid samples and include an 8-place sample block.

6.4 GC/MS Analytical System

The GC/MS analytical system includes a Gas chromatograph (Thermo Electron model Trace 1310 GC or equivalent) and a Mass Spectrometer (Thermo Electron model ISQ LT MS or equivalent). The system includes the following components:

6.4.1 Autosampler

The autosampler (Thermo Electron TriPlus model TP 100 or equivalent) is capable of delivering consistent volumes for the sample and the internal standard (IS). The autosampler can add IS prior to sample injection to the GC.

6.4.2 Gas Chromatograph (GC)

6.4.2.1 The GC oven is temperature programmable.

6.4.2.2 The sample inlet includes a large volume injector capable of injecting sample volume in the range of 1µL to 100µL, at a controlled rate (for example, 2µL/sec). The injector must be temperature-programmable and with the highest temperature ≥ 350 °C.

6.4.2.3 The column is a 30 meter long x 0.25mm ID fused silica capillary column coated with a cross-linked phenyl methyl silicone. Other column with a similar separation capability may also be used.

6.4.3 Mass Spectrometry (MS)

The MS is capable of producing electron impact spectra at +70 eV (nominal) and scanning the range of the specified quantification masses or m/z. A scan range of 40 to 500 m/z is adequate for PAH analysis.

6.4.3.1 Data Acquisition

The data acquisition is controlled by Xcalibur (Thermo Electron software), which searches any GC/MS runs for specific ions or reconstructed ions, and plots the intensity of the ions with respect to retention time or scan number.

7. Reagents and Materials

7.1 Gases

7.1.1 Compressed Helium, a carrier gas for GC/MS, at least 99.9995 % purity.

7.1.2 Compressed Nitrogen, at least 99.999% purity and free of PAHs contamination.

7.1.3 Compressed Zero Air, for ASE valve switching. It can be substituted by nitrogen gas specified in 7.1.2.

7.2 Reagents

7.2.1 Dichloromethane (DCM), at least 99.8% purity. The PAHs levels in this DCM must be less than the method detection limits and suitable for PAH trace analysis at or below the part-per-billion.

7.2.2 Acetone, at least 99.9% purity.

7.2.3 Hexane, at least 99% purity.

7.2.4 Methanol, at least 99% purity.

7.2.4 NIST Reference Materials. Any NIST reference materials containing the target analytes can be used for calibration, for example, SRM 1491.

7.2.5 Deionized water, at least 5Ωm resistance.

7.2.6 PAH standards, 10 µg/mL in Toluene, AccuStandard, Inc. Solutions from other sources are also acceptable as long as they contain PAHs of interest.

7.2.7 Deuterated PAHs

Deuterated PAHs are used as recovery and internal standards. The purity of these compounds shall be the highest whenever it is possible. The purities of these deuterated PAHs purchased from CDN are sufficient for the purposes.

7.3 Filters

Quartz fiber filter, glass fiber filter, and/or Teflon coated glass fiber filter such as Pallflex 7213 T60A20 are used for PM sampling.

7.4 PUF, polyurethane foam

7.5 Amberlite XAD4, surface area: 725m²/g and average pore diameter: 40A.

8. Preparation of Apparatus

8.1 ASE Extractor Cell Cleaning

- 8.1.1 Inspect the integrity of the extraction cell, such as the Teflon O-ring and seals. Replace them when necessary.
- 8.1.1 Clean the cells. Rinse the cell parts with DI water and acetone. Air dry the cell parts or oven dry at 100°C. Cell disposable filters can be cleaned with DCM via Soxhlet or by ASE extraction (dummy extraction) to ensure there is no PAH contamination from the disposable filters.
- 8.1.2 Sonicate the borosilicate glass balls (cell fillings) with DI water in a clean beaker for at least 30 minutes. Discard and drain the water. Repeat the procedure for the second time. Rinse the borosilicate glass balls with DI water, followed by acetone, hexane, and dichloromethane. Air dry or oven dry at 100°C in a PAH-free environment.

8.2 Glassware cleaning

All glassware can be cleaned by the following procedures:

- 8.2.1 Wash with detergent (Extran AP12, residue free or equivalent) and rinse with DI water during washing cycles in a washing machine.
- 8.2.2 Rinse the washed glassware three times with acetone, hexane and dichloromethane, in that order.
- 8.2.3 Dry the cleaned glassware in an oven at approximately 100 °C.
- 8.2.4 Alternative cleaning procedures can be used if the background concentrations for PAHs are below the method detection limits.

8.3 Soxhlet extractor cleaning

- 8.3.1 Assemble the Soxhlet extraction apparatus, as shown in Figure 2, with a thimble (if needed) in the extraction tube.
- 8.3.2 Set the temperature of the chiller at approximately 6 °C for condensing DCM.
- 8.3.3 Adjust the Variac to control the DCM distillation rate at approximately 3 cycles per hour. Continue distilling for at least 6 cycles.
- 8.3.4 Discard the solvent and dry the apparatus in the oven if desired.

9. Preparation of Sampling Materials

9.1 Filter cleaning

- 9.1.1 Filters for PAH analysis are pre-cleaned with DCM via Soxhlet extraction prior to sample collection.

- 9.1.2 Load the filters in the cleaned Soxhlet extraction tube and clean the filters according to the procedure in section 8.3. Continue distilling for 16 hours.
- 9.1.3 Remove the extraction tube. Purge the filters at room temperature with nitrogen gas till dry.
- 9.2 PUF cleaning
 - 9.2.1 PUF for PAH analysis is pre-cleaned with DCM via accelerating solvent extraction (ASE) following typical ASE conditions:
 - Pressure: 1500 psi
 - Temperature: 100°C
 - Static time: 5 min
 - Flush volume: 60%
 - Purge time: 100 sec
 - Static cycles: 3
 - Solvent: DCM
 - 9.2.2 Use vacuum oven to completely remove the residual solvent. The vacuum oven temperature is set at 31°C.
 - 9.2.3 If the size of the PUF is reduced too much, soak the PUF in DCM. Squeeze the PUF to remove excess solvent then use the vacuum oven to ensure complete removal of the DCM.
 - 9.2.4 This will allow the PUF to come back approximately to its original shape and size.
 - 9.2.5 Store the cleaned PUF in a jar at ~ 4°C.
- 9.3 XAD4 cleaning
 - 9.3.1 The XAD4 resin is cleaned with methanol and DCM by ASE following the conditions listed in section 9.2.1.
 - 9.3.2 Use vacuum oven to remove residual solvent at ~ 31°C.
 - 9.3.3 Store the XAD4 in a jar at ~ 4°C for later use.
 - 9.3.4 Analyze the last DCM wash to confirm the cleanness of the XAD4 resin.

10. Preparation of Standards for PAH Analysis

10.1 Preparation of Recovery Standards (RS)

The recovery standards (RS) consist of several deuterated PAHs for assessing the extraction efficiency. The target PAH species and their corresponding RS are listed in

Table 2. Other deuterated PAHs with similar structures may also be used depending upon their availability.

- 10.1.1 Accurately weigh each RS (approximately 0.01g).
- 10.1.2 Add all the weighed compounds to a 100 mL volumetric flask.
- 10.1.3 Dissolve the RS with DCM (< 100 mL) in the flask.
- 10.1.4 Add DCM to the 100mL mark to make a 100 µg/mL stock solution.
- 10.1.5 Pipet 1 mL of the 100 µg/mL stock solution, and dilute with DCM to the 100 mL mark. The dilution results in a 1 µg/mL solution and can be used for spiking the PUF/XAD/PUF samples before extraction.
- 10.1.6 Pipet 10 mL of the 1µg/mL solution and dilute to 100 mL to make a 100 ng/mL solution for further use in 10.2.
- 10.1.7 Pipet 4 mL of the 1µg/mL solution and dilute to 100 mL to make a 40 ng/mL recovery standard solution for spiking the filter sample prior to extraction.
- 10.1.8 Use the 40 ng/mL recovery standard to spike the filter samples if the sample extract needs to be concentrated more than four times. To calculate the volume of a recovery standard solution for spiking the filters, determine the desired concentrations of recovery standards in the extract including the final extract volume. The final concentrations should be within the linear range of the calibration standards. Use the 1 µg/mL recovery standard to spike the XAD/puff samples.
- 10.1.9 The integrity of the recovery standard requires frequent checking against calibration standards. The concentration of the recovery standard may change due to solvent evaporation, degradation, and any other reasons that can shorten its life.
- 10.1.10 Prior to extraction, check the concentration agreement between the deuterated PAH calibration standards and spiking solution by analyzing a freshly prepared QC solution of ~ 1 ng/mL from a spiking solution. If the results fall within 20 % of the expected values of each recovery standard, extract PM filters and XAD/PUF by following the procedures in section 11.
- 10.1.11 If the results exceed the 20% limit, correct the problems prior to extraction. If it is due to deterioration of deuterated compounds, prepare a new spiking solution and calibration standards from stock solution and repeat procedure listing in 10.1.10.

10.2 Preparation of PAH and deuterated PAH Calibration Standards

The calibration standards are obtained following a series dilution of the NIST Reference Material or other certified standard materials purchased from a qualified vendor. The calibration standards must contain at least the 23 certified PAHs. The PAH compounds along with their concentrations for each level of standard solutions are listed in Table 3. Other calibration standards including gravimetrically prepared standards from pure compounds are also acceptable.

- 10.2.1 The PAH calibration concentrations range from 0.05 to 5.0 ng/mL. The actual PAHs' concentrations may vary, depending on the accurate concentrations of the

standard materials being used. Table 3 lists the concentration levels used in this SOP.

- 10.2.2 Use a syringe to take out 500 μL of the PAH certified solution and transfer into to a 50 mL flask, then dilute with DCM to make a stock solution (100ng/mL).
- 10.2.3 Pipet 0.5, 1, 2, 4, and 5 mL of the stock solution and 0.5, 1, 2, 3 and 5 mL of the 100ng/mL recovery standard solution to 5 individual 100 mL flask, respectively, and dilute with DCM to the mark. The procedure will yield 5 calibration standards (Cal 2-6, respectively). Use Cal 2 solution to performing 1:10 dilution to make Cal 1 solution. The concentrations of these standards are listed in Table 3 for target PAHs and Table 4 for deuterated PAHs. The procedure and actual concentrations may vary due to availability of standard materials.

10.3 Preparation of deuterated PAH Internal standards (IS)

Five deuterated PAHs are selected as IS for the determination of PAH concentrations. Each analyte has one, in some cases two, of the IS assigned for quantification. The corresponding IS for each analyte is listed in Table 5.

- 10.3.1 The ISs used in this method are: acenaphthene-d10, phenanthrene-d10, pyrene-d10, chrysene-d12, and perylene-d12.
- 10.3.2 Should interference be found in the quantification ions of any IS (such as chrysene-d12 and pyrene-d10), use another deuterated PAH which has retention time close to that of the target analyte.
- 10.3.3 The following procedure is a typical way of preparing the IS solution. Weigh the four ISs individually (approximately 0.05 g) and add all weighed ISs into a 100 mL flask.
- 10.3.4 Dissolve the ISs with DCM. When the solid compound is dissolved completely, add more DCM to the 100 mL mark to make a $\sim 500 \mu\text{g/mL}$ stock solution.
- 10.3.5 Pipet 1 mL of the stock solution into a 100 mL flask and dilute with DCM to make a $5 \mu\text{g/mL}$ solution. Pipet 2 mL of the $5 \mu\text{g/mL}$ into a 100 mL flask and dilute with DCM to make a 100 ng/mL solution.
- 10.3.6 Pipet 5 mL of the 100 ng/ mL solution into a 100 mL flask and dilute with DCM to the mark to make an IS solution at a concentration of 5 ng/mL for each compound.

10.4 Preparation of QC standards

QC standards can be volumetrically diluted from a commercially available or a NIST SRM PAH standard. The concentrations of PAHs in the QC standards shall be in the mid range of the prepared GC/MS calibration standards. The suggested QC standard concentration is approximately 1.0 ng/mL.

- 10.4.1 A series of dilutions are needed. The actual dilution procedure is dependent upon the concentration of the solution. The following procedure is an example for diluting a 2 mg/mL PAH solution.
- 10.4.2 Use a syringe and take a 250 μL aliquot of the 2 mg/mL PAH standard into a 100 mL flask. Dilute the solution to the mark with DCM to make a $5.0 \mu\text{g/mL}$ solution.

- 10.4.3 Pipet 1mL of the 5 µg/mL solution and dilute to the 100 mL mark with DCM to make a 50ng/mL solution.
- 10.4.4 Continue pipetting 2 mL of the 50 ng/mL solution and dilute to the 100 mL mark with DCM to make a 1.0 ng/mL QC solution.

11. Procedures for PAH Quantification

11.1 Sample Extraction

The following section describes the procedures to operate the Dionex ASE350 and extract PAHs collected on PM filter or XAD/PUF media. The procedures should be modified if different sample collection media are used.

- 11.1.1 Visually examine the integrity of the cell parts, including Teflon O-ring, frits, and disposable filters. These parts provide seals needed for the pressurized and heated extraction cells.
- 11.1.2 Place the sample filter in the cell with tools (for example, forceps). Avoid direct contact between the forceps and the sample collection side of the filter to reduce potential losses of the analytes and to avoid cross contamination among samples. Should direct contact occur, rinse the forceps with acetone followed by hexane and DCM prior to next use, or use a clean forceps.
- 11.1.3 Spike the sample with the RS. The RS volume and concentration varies, depending upon whether the samples need to be concentrated prior to analysis. For example, spiking 100µL of a 1 µg/ mL RS solution to the XAD/PUF sample results in RS concentration of 1 ng/ mL in the 100 mL extract. Spiking 100 µL of 40 ng/mL RS solution to the filter sample results in RS concentration of 1.0 ng/mL in the 4 mL final extract solution if the volume of ASE extract, typically ~ 40 mL, is reduced 10 times.
- 11.1.4 Pack the cells with borosilicate glass balls. Do not over-pack the cells. Provide good seals for the cells to reduce potential extraction failure. Keep the threads and sealing surface clean.
- 11.1.5 Open the cylinder gas valves. Turn on the power of the ASE 350. Load the extract bottles. Use the ASE350 software to set up the extraction schedule and method (section 9.2.1). Start the extraction.
- 11.1.6 PAHs can be photosensitive and volatile. After extraction, measure the sample extract volume. Transfer the sample extracts to an amber bottle and store the bottle of extract in a refrigerator until analysis.
- 11.1.7 Extracting PAHs from the PUF/XAD/PUF can be performed directly if the Dionex extraction cell tube are used for packing the sampling media.
- 11.1.8 Complete the extraction by following the procedures list in 11.1.5-6.

11.2 Solvent Reduction

Solvent reduction is used to concentrate the extract after sample extraction in order to increase the method sensitivity. Estimate the final concentration prior to spiking the

recovery standards (RS). The expected RS final concentration shall be within the linear range of the calibration standards. Perform necessary dilution if the analytical results are outside the calibration range.

11.2.1 Nitrogen flow-down

The sample extract is concentrated by gently flowing nitrogen gas on top of the solution. This will evaporate the excess DCM and reduce the volume of the sample solution.

11.2.2 Connect the RapidVap N₂ dry evaporation system to the nitrogen gas source and allow nitrogen to flow through. The nitrogen can be either from a liquid nitrogen Dewar or from a zero nitrogen cylinder as long as the nitrogen gas is free of PAH contamination.

11.2.3 The following procedure is an example for concentrating a sample solution for ~10 times.

11.2.3.1 Weigh the weight of the dry glass tube equipped with the RapidVap N₂ dry evaporation system. Transfer ~20 mL of the sample solution to the glass tube and weigh the weight again. Calculate the original amount (mg) of solution that being transferred into the glass tube.

11.2.3.2 Start the nitrogen flow slowly to avoid splashing the solution. Set the temperature of the evaporation system at ~ 30 °C.

11.2.3.3 Stop nitrogen flow when the solution is below the 2.0 mL mark. Never allow the solution to dry out. Weigh the glass tube again and calculate the amount (mg) of solution remaining in the glass tube.

11.2.3.4 Calculate the exact concentration factor by determining the ratio between the amount (mg) of the original solution before the evaporation and the remaining solution after the evaporation. Other concentration factors are also acceptable as long as the PAH interference in the DCM is below the filter extract's detection limits.

11.2.3.5 The concentrated filter extract is ready for GC/MS analysis.

11.3 GC/MS Tuning

The following sections describe the procedures and conditions of conducting GC/MS analysis for PAHs.

11.3.1 Prior to the GC/MS analysis, the MS must pass its auto-tuning criteria with FC-43 (Perfluorotributylamine). The criteria include ion abundance, mass resolution, and mass calibration.

11.3.2 Auto-tune also automatically adjusts ion source parameters for optimum performance.

11.4 GC Operating Conditions

11.4.1 Typical autosampler parameters include:

Internal standard volume: 5.0 µL

Sample volume: 30.0 μL
Injection speed: 3.0 $\mu\text{L}/\text{sec}$
Syringe size: 100 μL

11.4.2 Typical GC operating conditions include the following:

GC column: 30 meter long x 0.25mm ID, DB-5 fused silica capillary column
Oven method: initial temperature: 40 $^{\circ}\text{C}$ for 2 minutes
Temperature ramping rate: 10 $^{\circ}\text{C}/\text{min}$
Final temperature: 300 $^{\circ}\text{C}$ and hold for 10 mins.
Carrier gas: helium with a constant flow rate at 1.20 mL/min.

11.4.3 Sample introduction is operated on a Programmable temperature and volume (PTV) large volume injection mode. A typical sample inlet method includes an initial condition and four phases. The Injection phase allows the sample to be introduced into the inlet at a low temperature with a vent valve open to reduce the excess DCM vapor volume and concentrates the analytes inside the injector. Evaporation phase allows further reduction of the solvent volume. Transfer phase transfers the sample to the GC column at a higher temperature. Finally, the inlet undergoes the Cleaning phase. The following is a typical sample injection method:

Initial conditions:

Base temperature: 35 $^{\circ}\text{C}$
Split flow rate: 20 mL/min
Split Ratio: 6.7.
Splitless time: 0.5 min

Injection phase:

Injection time: 0.65 min
Vent flow: 20 mL/min

Evaporation phase:

Evaporation temperature ramping rate: 14.5 $^{\circ}\text{C}/\text{sec}$
Evaporation temperature: 50 $^{\circ}\text{C}$
Evaporation time: 0.15 min

Transfer phase:

Transfer temperature ramping rate: 2.5 $^{\circ}\text{C}/\text{sec}$
Transfer temperature: 300 $^{\circ}\text{C}$
Transfer time: 0.5 min

Cleaning Phase:

Clean rate: 14.5 $^{\circ}\text{C}/\text{min}$
Clean temperature: 325 $^{\circ}\text{C}$
Clean time: 25.0 min
Clean flow: 80 mL/min

11.5 MS Operation Conditions

11.5.1 The mass spectrometer is operated on a positive electron impact (EI+) mode. The typical settings for the MS include:

Source temperature: 250 °C
GC interface temperature: 300 °C
Emission current 50 µA
Electron Energy: 70.0 V (nominal)

11.5.2 The MS can be operated under either Selective Ion Monitoring Mode (SIM) or Full Scan Mode. SIM is usually used for quantification and Full Scan can be used for both compound identification and quantification. The primary (or quantifier) and secondary (or qualifier) ions for quantifying each analyte and each internal standard in the SIM mode are listed in Tables 5 and 6, respectively. For the Full Scan, the mass scan range is 50 to 500 amu at a scan rate of one second per scan.

11.6 PAH Analysis

11.6.1 Prior to PAH analysis, allow the calibration standards (or RS if assessing recovery efficiency), the extract samples, the IS, and the QC solutions to reach room temperature.

11.6.2 Clean and fill the washing vials with DCM. All the vial septa for GC analysis must be Teflon to avoid PAH contamination.

11.6.3 Tune the mass spectrometer according to manufacturer's instructions, as described in Section 11.3.

11.6.4 Set GC/MS operation conditions in accordance with sections 11.4 and 11.5.

11.6.5 Load the DCM blank, calibration standards, sample extracts, IS, and QC samples to the GC/MS autosampler. Start the analysis with the DCM blank followed by calibration standards to evaluate GC/MS background and performance.

11.6.6 Analyze at least one QC and one duplicate sample for every 10 samples.

11.6.7 Set up the sequence for sample analysis.

11.7 Calibration Curve Construction

Construct the calibration curves using SIM and internal standard method according to Section 13.1 and 13.2. Check the r^2 of linear regression. If the r^2 is less than the specified value (0.99), recalibrate the GC/MS system or correct any problem.

12. Quality Control (QC)

12.1 Calibration

12.1.1 GC/MS requires frequent calibration. Calibration is needed when starting an analysis or any QC failure, as described in section 12.4.

- 12.1.2 Check the correlation r^2 value for each PAH calibration curve. The r^2 value must be at least 0.99.
- 12.1.3 Failure to achieve the linearity requirement will require re-calibrating the GC/MS. If it fails after recalibration, the problems causing the failures should be corrected before any analysis can be reported.

12.2 Method detection limit and reporting limit

Method detection limit (MDL) for a sample extract is defined as follows:

$$MDL(i) = t * SD (i)$$

where:

$MDL(i)$ = method detection limit of PAH i

t = student's t value associated with a 98% confidence interval (section 12.3.1)

$SD (i)$ = standard deviation of a replicate analysis of the lowest concentration standard for PAH component i

- 12.2.1 The Student's t value is dependent upon the degrees of freedom associated with the analysis. The degree of freedom of the analysis is equal to the number of replicate measurements, n , of the lowest concentration standard, minus one. An abbreviated table of values of t associated with a 98% confidence interval is shown below^(ref:4):

Degrees of Freedom (n-1)	t-value
4	3.7
5	3.4
6	3.1
7	3.0

- 12.2.2 The MDL for PAHs should be verified annually or when the method is modified with potential influence on the detection limits. It is carried out by repeating the lowest concentration of PAH standards analysis.
- 12.2.3 The limit of each PAH's reporting limit (RL) for a sample extract is set between 0.03 and 0.1 ng/mL. However, Naphthalene and methyl-naphthalenes are likely subjected to background interferences, their MDLs need to be examined carefully.
- 12.2.4 If the calculated MDL for sample extract is greater than the RL and the causes of the elevated MDLs are due to environmental influences, the reporting limit, may be reset as the calculated MDL.
- 12.2.5 The typical reporting limit for PAH mass on the PM filter is in the range of 0.12 ~ 0.40 ng per filter based on a final extraction volume of 4 mL. The reporting limits may vary due to a variation of the final extract volume and analytical conditions.
- 12.2.6 The typical reporting limit for PAH mass on the XAD/PUF cartridge is in the range of 3.0 ~ 10.0 ng per sample and may vary due to a variation of the final extract

volume and analytical conditions.

12.3 QC sample analysis

A QC sample is analyzed after calibration, after every 10 samples, and at the end of a sample set. If a QC failure occurs, the sources of the problems must be corrected before conducting any additional sample analysis.

12.3.1 Not all PAHs are monitored for quality control purposes. The choice of QC compounds depends upon the composition of commercially available standards and the retention times of the PAHs. This SOP uses the following PAHs to determine the performance of the GC/MS: naphthalene, acenaphthylene, fluoranthene, benzo[b]fluoranthene, benzo[a]pyrene, and benzo[g,h,i]perylene.

12.3.2 Analyze the QC standards to determine the concentrations of these QC compounds.

12.3.3 The QC limit is 20% from the calculated concentration for each component in the QC standards. The QC is considered "pass" if the QC result is within 20% of the calculated concentration. Any result outside 20% is considered as a QC "failure". Problems causing a QC failure must be determined and corrected before further samples are analyzed.

12.4 Blank sample Analysis

12.4.1 A DCM sample must be analyzed at the beginning of the sequence of the analysis. The PAH level in the DCM must be below the PAH's reporting limit before proceeding to instrumental calibration or sample analysis.

12.4.2 Three trip blank samples must be analyzed in each sequence of the analysis. The PAH level in the trip blank samples must not exceed the PAH levels determined in the field blank samples or any field samples.

12.4.3 The PAH levels in the trip blank samples must be subtracted from the field blank and field samples if the PAH levels in the trip blank samples are higher than the reporting limit.

12.5 Replicate Analysis

For each set of samples, perform a replicate sample analysis at least once or after every 10 sample analyses. The difference of the two replicate results for QC compounds listed in section 12.3.1 must be within 10% of their mean.

12.6 PAH identification

Identify the various PAHs in Table 1 from their retention times and mass spectra. Confirmation of a PAH can be obtained by matching the retention time and full scan MS of the GC/MS analysis with a calibration analysis. When comparing the criteria below must be met:

12.6.1 All characteristic ion peaks for the same PAH must have retention times within 0.02 min of one another.

12.6.2 The area ratio of the secondary to primary ion for the PAH must be within 5% when compared to the same ion area ratio from a calibration standard at approximately the same concentration. If such agreement cannot be made,

choose secondary ions for calibration assuming the secondary ions are free of interferences or flag for "possible interference" in the final reports.

- 12.6.3 The GC retention time (from total ion chromatogram) must be within 0.03 min of that obtained for the same calibration standard.

13. Calculation

13.1 Calibration:

13.1.1 After analysis of the standards, integrate the characteristic ion's (Table 5 and 6) peak area for each calibration compound and internal standard at the expected retention time. Figure 3 shows the chemical structure and the mass spectrum of each PAH compound in this analysis (Table 1). Figure 4 shows the extracted ion chromatogram (EIC) of the calibration standard PAHs.

13.1.2 Plot the response area ratio $R_{sp}(i)$ (y-axis) versus the amount ratio $Amt(i)$ (x-axis) to generate calibration curves for each PAH compound listed in Table 1:

Where:

$$R_{sp}(i) = A(i)/A(is)$$

$A(i)$ = area of compound i , and
 $A(is)$ = area of internal standard.

and $Amt(i) = C(i)/C(is)$

$C(i)$ = concentration of compound i in the calibration standard,
 $C(is)$ = concentration of the IS.

13.2 Linear least squares fit calculation

Obtain the linear least squares fit from the calibration curves. For each PAH i , the linear least squares fit can be expressed in the following form:

$$R_{sp}(i) = m(i) \times Amt(i) + b(i)$$

where

$m(i)$ = slope of linear equation for PAH i

$b(i)$ = intercept at the y-axis

13.3 PAH concentration calculation

From the equations in Sections 13.1.2 and 13.2, calculate the concentration of each PAH, $C(i)$, in ng/mL, in the filter extract using the area ratio ($R_{sp}(i)$) of the area for the sample of the PAH to that of the IS.

$$C(i) = [(A(i)/A(is) - b(i)) / m(i)] \times C(is)$$

13.4 Recovery efficiency calculation

The recovery efficiency can be calculated by comparing the original spiking deuterated PAH concentrations to the concentrations obtained from a GC/MS.

13.4.1 Follow the procedures described in sections 13.1-13.3 to determine the concentration of each deuterated PAH in the filter extract.

13.4.2 Calculate the original theoretical concentration of each deuterated PAH in the filter extract as follows:

$$Theo(i) = RSpike(i) \times V(inj) / V(fext)$$

where

Theo(i) = Calculated original theoretical concentration of deuterated PAH component *i* in the filter extract,

RSpike(i) = Concentration of deuterated PAH component *i* in spike solution,

V(inj) = volume of spike solution injected on the filter prior to extraction, and

V(fext) = final volume of filter extract

13.4.3 Calculate the recovery efficiency for each deuterated PAH as follows:

$$\% Recovery(i) = [C(i) / Theo(i)] \times 100\%$$

where

% Recovery(i) = recovery percentage of deuterated PAH component *i*

C(i) = concentration of deuterated PAH component *i* obtained from the GC/MS, and

Theo(i) = calculated original theoretical concentration of deuterated PAH component *i* in the filter extract.

13.4.4 If the calculated efficiency is less than 50% or higher than 150%, a remark or flag "recovery low" or "recovery high" must be noted on the analysis report.

13.5 Calculate the PAH mass in the sample

13.5.1 Calculate individual PAH mass

The individual PAH mass measured from each sample can be calculated by multiplying the extract volume and corrected for its extraction recovery efficiency.

13.5.2 The PAH mass can be calculated as follows:

$$PAH(i) = C(i) \times V(fext) / \% Recovery(i)$$

where

PAH(i) = measured mass (in ng) of PAH component *i* on the filter or XAD/PUF,

C(i) = concentration(in ng/mL, or ppb[w/v]) of PAH *i* obtained from GC/MS analysis,

V(fext) = Final sample's extraction volume (in mL) prior to GC/MS analysis, and

% Recovery(i) = recovery percentage of PAH *i*.

13.5.3 For PAHs that have recovery efficiency greater than 100%, the recovery efficiency of 100% will be used to calculate the PAH mass.

14. Report

The report for each filter sample shall contain the identified PAH names and their individual measured masses, flags for high or low recovery efficiency, and flags for any possible interference for concentration determination. An example of the report is shown in Figure 5.

Any results below the reporting limit will be reported as "<RL".

15. Reference Documents

1. California Environmental Protection Agency, Air Resources Board, Method 429: Determination of Polycyclic Aromatic Hydrocarbon (PAH) Emissions from Stationary Sources.
2. US Environmental Protection Agency, Analytical Method for the Analysis of Semi-volatile Organic Compounds, Exhibit D
3. Desert Research Institute Compendium Method TO-13A: Determination of Polycyclic Aromatic Hydrocarbons (PAHs) in Ambient Air Using gas Chromatography/mass Spectrometry (GC/MS).
4. Harris, Daniel C., "Quantitative Chemical Analysis", *W.H. Freeman & Co.*, 4th ed., 1995.

16. Revision History

1. SOP NO. MLD 144 V 1.0 was approved in November 2006.
2. To reflect the CARB division realignments, SOP NO. MLD 144 changes its convention to MV-AEROSL-144. The version 2.0 was approved in June, 2014.
 - a. The procedure for determination the semi-volatile PAHs were added to Version 2.0.
 - b. Version 2.0 also incorporates editorial changes.
3. MV-AEROSL-144 V 2.1 was approved in June, 2019 to incorporate editorial changes and method modifications due to the changes of GC-MS system and sample concentrator.
4. MV-AEROSL-144 V 2.2 was approved in September, 2019 to incorporate changes of the Division name.

Figure 1. Flowchart of Filter Sample Analysis




Figure 2. Soxhlet Extraction Apparatus










Figure 3. The molecular structure and the mass spectrum of the PAH compounds analyzed by GC/MS in this method.







Acenaphthene









2,3,5-trimethylnaphthalene









Fluorene
















Figure 4. The Extracted Ion Chromatogram (EIC) of the Calibration Standard PAHs




Figure 5. Sample report of the PAH analysis using GC/MS
PAH CONTENTS IN PM SAMPLE

Chemical Analysis and Emissions Research Branch
 Emission Compliance, Automotive Regulations and Science Division
 California Air Resources Board

Sample ID: TO2208B
 Data File: 052406 TO2208B
 Curr Data Path: C:\Xcalibur\Data\
 Operator: XAD
 Acquisition Date: 05/25/06 05:13:43
 Inst Method: C:\Xcalibur\methods\SIM_SRM1491.meth

Target Analyte	Calculated Amount (ng/sample)	Remarks
010_Naphthalene	7990.913	
015_2-methylnaphthalene	1085.267	
018_1-methylnaphthalene	810.228	
050_2,6-dimethylnaphthalene	50.105	
070_Acenaphthylene	792.821	
090_Acenaphthene	112.318	
100_2,3,5-trimethylnaphthalene	8.368	
110_Fluorene	383.010	
210_Phenanthrene	264.929	
220_Anthracene	115.214	
230_1-methylphenanthrene	20.618	
240_Fluoranthene	68.045	
310_Pyrene	160.794	
320_Benz[a]anthracene	<RL	
410_Chrysene	<RL	
420_Benzo[b]fluoranthene	<RL	
430_Benzo[k]fluoranthene	<RL	
440_Benzo[e]pyrene	<RL	
450_Benzo[a]pyrene	<RL	
470_Perylene	<RL	
480_Indeno[1,2,3-cd]pyrene	<RL	
490_Dibenz[a,h]anthracene	<RL	
495_Benzo[ghi]perylene	<RL	

Table 1. Target PAHs and formulae

	Compound	Formula
1	Naphthalene	C ₁₀ H ₈
2	2-Methylnaphthalene	C ₁₁ H ₁₀
3	1-Methylnaphthalene	C ₁₁ H ₁₀
4	2,6-Dimethylnaphthalene	C ₁₂ H ₁₂
5	Acenaphthylene	C ₁₂ H ₈
6	Acenaphthene	C ₁₂ H ₁₀
7	2,3,5-Trimethylnaphthalene	C ₁₃ H ₁₄
8	Fluorene	C ₁₃ H ₁₀
9	Phenanthrene	C ₁₄ H ₁₀
10	Anthracene	C ₁₄ H ₁₀
11	1-Methylphenanthrene	C ₁₅ H ₁₂
12	Fluoranthene	C ₁₆ H ₁₀
13	Pyrene	C ₁₆ H ₁₀
14	Benz[a]anthracene	C ₁₈ H ₁₂
15	Chrysene	C ₁₈ H ₁₂
16	Benzo[b]fluoranthene	C ₂₀ H ₁₂
17	Benzo[k]fluoranthene	C ₂₀ H ₁₂
18	Benzo[e]pyrene	C ₂₀ H ₁₂
19	Benzo[a]pyrene	C ₂₀ H ₁₂
20	Perylene	C ₂₀ H ₁₂
21	Indeno[1,2,3-cd]pyrene	C ₂₂ H ₁₂
22	Dibenz[a,h]anthracene	C ₂₂ H ₁₄
23	Benzo[ghi]perylene	C ₂₂ H ₁₂

Table 2. Target PAHs and their Recovery Standard

Compound	Recovery standard
Naphthalene	Naphthalene - d8
2-Methylnaphthalene	2-Methylnaphthalene-d10
1-Methylnaphthalene	2-Methylnaphthalene-d10
2,6-Dimethylnaphthalene	2,6-Dimethylnaphthalene-d12
Acenaphthylene	Acenaphthylene-d8
Acenaphthene	Acenaphthylene-d8
2,3,5-Trimethylnaphthalene	Acenaphthylene-d8
Fluorene	Acenaphthylene-d8
Phenanthrene	Anthracene- d10
Anthracene	Anthracene- d10
1-Methylphenanthrene	Anthracene- d10
Fluoranthene	Fluoranthene-d10
Pyrene	Fluoranthene-d10
Benz[a]anthracene	Benz[a]anthracene-d12
Chrysene	Benz[a]anthracene-d12
Benzo[b]fluoranthene	Benzo[b]fluoranthene-d12
Benzo[k]fluoranthene	Benzo[k]fluoranthene-d12
Benzo[e]pyrene	Benzo[a]pyrene-d12
Benzo[a]pyrene	Benzo[a]pyrene-d12
Perylene	Benzo[a]pyrene-d12
Indeno[1,2,3-cd]pyrene	Indeno[1,2,3-cd]pyrene-d12
Debenz[a,h]anthracene	Debenz[a,h]anthracene-d14
Benzo[ghi]perylene	Benzo[ghi]perylene-d12

Table 3. PAH Calibration standards

Unit: ng/mL

Compound	Cal 1	Cal 2	Cal 3	Cal 4	Cal 5	Cal 6
Naphthalene	0.050	0.499	0.998	1.996	3.992	4.990
2-Methylnaphthalene	0.049	0.495	0.989	1.978	3.956	4.945
1-Methylnaphthalene	0.050	0.499	0.999	1.998	3.996	4.995
2,6-Dimethylnaphthalene	0.050	0.499	0.999	1.998	3.996	4.995
Acenaphthylene	0.050	0.500	1.000	2.000	4.000	5.000
Acenaphthene	0.050	0.503	1.006	2.012	4.024	5.030
2,3,5-Trimethylnaphthalene	0.050	0.503	1.005	2.010	4.020	5.025
Fluorene	0.050	0.504	1.007	2.014	4.028	5.035
Phenanthrene	0.050	0.496	0.991	1.982	3.964	4.955
Anthracene	0.050	0.504	1.009	2.018	4.036	5.045
1-Methylphenanthrene	0.050	0.5	0.999	1.998	3.996	4.995
Fluoranthene	0.050	0.503	1.005	2.010	4.020	5.025
Pyrene	0.050	0.502	1.003	2.006	4.012	5.015
Benz[a]anthracene	0.050	0.501	1.002	2.004	4.008	5.010
Chrysene	0.050	0.498	0.996	1.992	3.984	4.980
Benzo[b]fluoranthene	0.050	0.503	1.006	2.012	4.024	5.030
Benzo[k]fluoranthene	0.050	0.502	1.003	2.006	4.012	5.015
Benzo[e]pyrene	0.050	0.502	1.004	2.008	4.016	5.020
Benzo[a]pyrene	0.050	0.498	0.996	1.992	3.984	4.980
Perylene	0.050	0.502	1.003	2.006	4.012	5.015
Indeno[1,2,3-cd]pyrene	0.050	0.501	1.002	2.004	4.008	5.010
Debenz[a,h]anthracene	0.050	0.499	0.998	1.996	3.992	4.990
Benzo[ghi]perylene	0.049	0.490	0.980	1.960	3.920	4.900

Table 4. PAH Recovery Standard Calibration standards

Unit: ng/mL

Compound	Cal 1	Cal 2	Cal 3	Cal 4	Cal 5	Cal 6
Naphthalene-d8	0.050	0.50	1.00	2.00	3.00	5.00
2-Methylnaphthalene-d10	0.050	0.50	1.00	2.00	3.00	5.00
2,6-Dimethylnaphthalene-d12	0.050	0.50	1.00	2.00	3.00	5.00
Acenaphthylene-d8	0.050	0.50	1.00	2.00	3.00	5.00
Anthracene-d10	0.050	0.50	1.00	2.00	3.00	5.00
1-Methylphenanthrene	0.050	0.50	1.00	2.00	3.00	5.00
Fluoranthene-d10	0.050	0.50	1.00	2.00	3.00	5.00
Benz[a]anthracene-d12	0.050	0.50	1.00	2.00	3.00	5.00
Benzo[b]fluoranthene-d12	0.050	0.50	1.00	2.00	3.00	5.00
Benzo[k]fluoranthene-d12	0.050	0.50	1.00	2.00	3.00	5.00
Benzo[a]pyrene-d12	0.050	0.50	1.00	2.00	3.00	5.00
Indeno[1,2,3-cd]pyrene-d12	0.050	0.50	1.00	2.00	3.00	5.00
Debenz[a,h]anthracene-d14	0.050	0.50	1.00	2.00	3.00	5.00
Benzo[ghi]perylene-d12	0.050	0.50	1.00	2.00	3.00	5.00

Table 5. Target PAHs and their Internal Standard (IS)

Target PAHs	Internal standard
Naphthalene	Acenaphthene-d10
2-Methylnaphthalene	Acenaphthene-d10
1-Methylnaphthalene	Acenaphthene-d10
2,6-Dimethylnaphthalene	Acenaphthene-d10
Acenaphthylene	Acenaphthene-d10
Acenaphthene	Acenaphthene-d10
2,3,5-Trimethylnaphthalene	Acenaphthene-d10
Fluorene	Acenaphthene-d10
Phenanthrene	Phenanthrene-d10
Anthracene	Phenanthrene-d10
1-Methylphenanthrene	Phenanthrene-d10
Fluoranthene	Pyrene-d10
Pyrene	Pyrene-d10
Benz[a]anthracene	Pyrene-d10
Chrysene	Chrysene-d12 or Perylene-d12
Benzo[b]fluoranthene	Chrysene-d12 or Perylene-d12
Benzo[k]fluoranthene	Chrysene-d12 or Perylene-d12
Benzo[e]pyrene	Chrysene-d12 or Perylene-d12
Benzo[a]pyrene	Chrysene-d12 or Perylene-d12
Perylene	Chrysene-d12 or Perylene-d12
Indeno[1,2,3-cd]pyrene	Chrysene-d12 or Perylene-d12
Debenz[a,h]anthracene	Chrysene-d12 or Perylene-d12
Benzo[ghi]perylene	Chrysene-d12 or Perylene-d12

Table 6. Target PAHs and their Characteristic Ions for Quantification

Target PAHs	Primary Ion	Secondary ion(s)
Naphthalene	128	127
2-Methylnaphthalene	142	141
1-Methylnaphthalene	142	141
2,6-Dimethylnaphthalene	156	155,154
Acenaphthylene	152	151
Acenaphthene	153	154,152
2,3,5-Trimethylnaphthalene	155	170,165
Fluorene	166	165
Phenanthrene	178	177,176
Anthracene	178	177,176
1-Methylphenanthrene	192	191
Fluoranthene	202	201,200
Pyrene	202	201
Benz[a]anthracene	228	229
Chrysene	228	229
Benzo[b]fluoranthene	252	250,253
Benzo[k]fluoranthene	252	250,253
Benzo[e]pyrene	252	250,253
Benzo[a]pyrene	252	250,253
Perylene	252	250,253
Indeno[1,2,3-cd]pyrene	276	277
Debenz[a,h]anthracene	278	276
Benzo[ghi]perylene	276	277

Table 7. Characteristic Ions for Internal and Recovery Standards

Compound	Type*	Primary Ion	Secondary ion(s)
Naphthalene-d8	RS	136	135
2-Methylnaphthalene-d10	RS	152	150
2,6-Dimethylnaphthalene-d12	RS	168	150
Acenaphthylene-d8	RS	160	159
Acenaphthene-d10	IS	164	163
Phenanthrene-d10	IS	188	189
Anthracene-d10	RS	188	189
Fluoranthene-d10	RS	212	211
Pyrene-d10	IS	212	211
Benz[a]anthracene-d12	RS	240	241
Chrysene-d12	IS	240	241
Benzo[b]fluoranthene-d12	RS	264	265
Benzo[k]fluoranthene-d12	RS	264	265
Benzo[a]pyrene-d12	RS	264	265
Perylene-d12	IS	264	265
Indeno[1,2,3-cd]pyrene-d12	RS	288	289
Debenz[a,h]anthracene-d14	RS	292	288
Benzo[ghi]perylene-d12	RS	288	289

RS: Recovery Standard**IS: Internal Standard**