

M435 Implementation Guidance Document Topics

- Introduction
- Applicability
- Sampling Practices
- Sample Processing Practices
- Laboratory Sample Analysis
- Quality Control

2

Introduction Background

- 1990 ARB Asbestos Airborne Toxic Control Measure (ATCM) for Surfacing Applications.
- 2001 ARB Asbestos ATCM for Construction, Grading, Quarrying, and Surface Mining Operations.
- 2007 ARB M435 Interlaboratory Study (ILS).
 - Sample processing/analytical procedures varied.
 - Can affect reported asbestos content.

3

Introduction Purpose

- Assist asbestos stakeholders in the application and performance of M435.
- Clarify procedures and recommend good field sampling and laboratory practices.
- Result in more accurate and repeatable M435 asbestos content measurements.


4

Applicability of M435

- Surfacing ATCM - M435 is the **required test method** to determine the asbestos content of surfacing aggregate materials.
- Construction ATCM - M435 is referenced as a **laboratory test method** to determine the asbestos content of bulk samples.
- M435 random sampling plan is not applicable for asbestos ATCM exemption requirements.

5

Sampling Practices for Aggregate Materials Sources and Sampling Design



6

Sampling Practices for Aggregate Materials Equipment and Procedures

- Storage piles - sloughing effect.
 - Sampling tubes, round point shovels, front loaders.
- Conveyor belts – less sloughing.
 - Automated or manual sampling using templates.
- Aggregate-covered surfaces.
 - Manual or automated augers, shovel, or other suitable equipment.



7


Sampling Practices for Aggregate Materials Sample Documentation



- Complete sample description per M435.
- Sample acceptance criteria.
- Chain of custody documentation.
- Laboratory information management system.

8

Laboratory Sample Processing Practices M435 Sample Processing



9

Laboratory Sample Processing Practices Drying

- Purpose of drying: to remove moisture that would hinder complete pulverization.
- Recommend standardization of laboratory drying procedures (e.g., temperature, time).



10

Laboratory Sample Processing Practices Crushing

- Sample must be crushed to a nominal size of <3/8 inch (~0.95 centimeter).
- Recommend use of jaw crushers:
 - Reliable.
 - Timely.
- Not recommended:
 - Use of hammers.
 - Restricting rock sample sizes.
 - Discarding large rock samples.



11

Laboratory Sample Processing Practices Sequence of Post-crushing Procedures

- Post-crushing processing procedures:
 - a) Sample size reduction.
 - b) Homogenization.
 - c) Pulverization.
- The sequence of these procedures affects the representativeness of the powder analyzed.
- The sequence also depends on what equipment are present in the laboratory.

12

Laboratory Sample Processing Practices Sequences of Post-crushing Procedures

Available Post-crushing Equipment	Step 1	Step 2	Step 3
Braun Mill + Mixer (Sequence 1)	Pulverize entire crushed sample.	Use mixer to homogenize entire powdered sample.	Manually take 1 pint for analysis.
Braun Mill (no Mixer) (Sequence 2)	Pulverize entire crushed sample.	Manually homogenize powdered sample.	Manually take 1 pint for analysis.
Shatterbox (SB), Ball Mill (BM), or Freezer Mill (FM) + Mixer (Sequence 3)	Use mixer to homogenize entire crushed sample.	Manually take 1 pint crushed sample for pulverization.	Pulverize 1 pint crushed sample for analysis.
SB, BM, or FM (no Mixer) (Sequence 4)	Riffle split entire crushed sample and take 1 pint aliquot.	Pulverize 1 pint crushed sample.	Manually homogenize powdered sample for analysis.
M435: Braun Mill or Equivalent (no Mixer)	Riffle split entire crushed sample and take 1 pint aliquot.	Pulverize 1 pint crushed sample for analysis.	

Test Method 435

13

Laboratory Sample Processing Practices Sequences of Post-crushing Procedures

Recommended

Available Post-crushing Equipment	Step 1	Step 2	Step 3
Braun Mill + Mixer (Sequence 1)	Pulverize entire crushed sample.	Use mixer to homogenize entire powdered sample.	Manually take 1 pint for analysis.
Braun Mill (no Mixer) (Sequence 2)	Pulverize entire crushed sample.	Manually homogenize powdered sample.	Manually take 1 pint for analysis.
Shatterbox (SB), Ball Mill (BM), or Freezer Mill (FM) + Mixer (Sequence 3)	Use mixer to homogenize entire crushed sample.	Manually take 1 pint crushed sample for pulverization.	Pulverize 1 pint crushed sample for analysis.
SB, BM, or FM (no Mixer) (Sequence 4)	Riffle split entire crushed sample and take 1 pint aliquot.	Pulverize 1 pint crushed sample.	Manually homogenize powdered sample for analysis.
M435: Braun Mill or Equivalent (no Mixer)	Riffle split entire crushed sample and take 1 pint aliquot.	Pulverize 1 pint crushed sample for analysis.	

Test Method 435

14

Laboratory Sample Processing Practices Pulverization

- Recommend the use of the Braun Mill for pulverization.
- ILS results - Braun Mill produced powder with:
 - No leftover chunks.
 - Not over-pulverized.
 - Majority of particles: <math><75\mu\text{m}</math> diameter.



15

Laboratory Sampling Processing Practices Pulverization

- Equivalency of other pulverizers to the Braun Mill.
- M435-specific pulverization protocol for this equipment.
- Acceptable particle size distribution and equivalent size characteristics.



16

Laboratory Sample Processing Practices Pulverization

The recommended PSD:

- $\geq 98\%$ of the powder passes through the 250- μm mesh sieve.
- 75- to 250- μm size fraction is 40 to 50% of the total mass of sample processed.
- $< 75\text{-}\mu\text{m}$ size fraction is 50 to 60% of the initial sample mass.



17

Laboratory Sample Processing Practices Homogenization

- Increases powder homogeneity.
- Increases the likelihood that a representative aliquot of the field sample is analyzed.
- Increases accuracy and repeatability.
- Recommend:
 - use of a 3-dimensional mixer.



18

Laboratory Sample Processing Practices Sample Size Reduction of Crushed Sample


- Even number of equal width chutes.
- At least 8 chutes for coarse aggregate, or 12 chutes for fine aggregate.
- Minimum width of the individual chutes should be about 50% larger than the largest particles in the sample.



19

Laboratory Sample Analysis Procedures Principles of M435 Asbestos Identification

- Morphology requirements.
- Optical characteristics determination by polarized light microscopy (M435 Table 3).



20

Laboratory Sample Analysis Procedures PLM Limits of Resolution

- Fine particles and fibers $\leq 2\mu\text{m}$ in length and $\leq 0.15\mu\text{m}$ in thickness.
- Smaller asbestos fibers may still be present.
- Other analytical techniques are not part of M435:
 - X-ray diffractometry (XRD).
 - Scanning electron microscopy (SEM).
 - Transmission Electron Microscopy (TEM).
- The Department of Toxic Substances Control uses a tiered analytical approach of PLM followed by TEM.

21


Laboratory Sample Analysis Procedures Asbestos Quantification

- Required testing volume:
 - 1 pint powdered aliquot.
- Recommended powder mass:
 - 5 mg per PLM slide, particle loading of $\sim 30\%$.
- Fiber identification requirements:
 - Length-to-width aspect ratio of $\geq 3:1$ **and**
 - Asbestos optical properties (M435 Table 3).
 - Only asbestos characteristics, as described in M435 Table 3, should be used for asbestos identification.

22

Laboratory Sample Analysis Procedures Asbestos Quantification

- Recommended point-count reticle:
 - Standard crosshair reticle.
- Additional objective lens recommended:
 - 20X PLM objective to verify optical properties.
- Increase in points counted:
 - Multiples of 400 (e.g., 800, 1,200, 1,600, etc.).
 - Additional points counted may increase false negative errors.



23

Laboratory Sample Analysis Procedures Exceptions for Point-counting

- Exception I – sample contains no asbestos.
 - No asbestos fibers in 10 fields of view (FOVs) in each of 3 PLM slides (30 FOVs total).
 - Report no asbestos was found by visual technique.
- Exception II – sample contains $>10\%$ asbestos.
 - Asbestos content exceeds 10% from 3 PLM slides.
 - Report asbestos content exceeds 10% by visual technique.

24

Quality Control Sampling

- Sampling equipment cleanliness.
- Equipment cleaning protocol.
- Integrity of field samples.
- Protection of sample identity.
- Field log of M435 samples, including sampling details.



25

Quality Control Processing

- Sample chain of custody.
- Laboratory information management system.
- Written SOP specific for M435.
- Processing blanks alongside regular field samples.
- Calibration of processing equipment.
- Particle size calibration check.



26

Quality Control Analysis

- Microscope alignment.
- Refractive index liquid calibration.
- Asbestos proficiency training.
- Blind analytical replicates.
- Instrument cross checks.
- M435 laboratory protocol validation.
- Performance evaluation slides.
- Documentation of results.



27

California Environmental Protection Agency
Air Resources Board

For questions and comments, please contact:

Jeff Wright
Manager, DASPS
jeffrey.wright@arb.ca.gov
916.322.7055

Rebecca D. Neumann, Ph.D., P.G.
Air Pollution Specialist
rebecca.neumann@arb.ca.gov
916.324.1145

Visit the Test Method 435 website and join the Asbestos List Serve:
<http://www.arb.ca.gov/toxics/asbestos/tm435/tm435.htm>

ARB Monitoring and Laboratory Division
Data Analysis and Special Projects Section (DASPS)
P.O. Box 2815, Sacramento CA 95812

28